


## Analysis Computer (BMB)

### INSTRUMENT SPECIFICATION SHEET



<b>Location</b>	Room 134, Bio21 Institute	
<b>Booking</b>	<a href="http://www.bio21.org/booking/month.php?user=IMAGEANAL&amp;rid=64">http://www.bio21.org/booking/month.php?user=IMAGEANAL&amp;rid=64</a>	
<b>Base</b>	Dell Precision T7600	
<b>Operating System</b>	Windows 7 Enterprise SP1 64 bit	
<b>Memory</b>	32GB	
<b>Monitor</b>	Dell U2312HM 23"	
<b>Video Card</b>	2.5 GB NVIDIA Quadro 5000	
<b>Hard Drive</b>	512GB	
<b>Processor</b>	Dual Intel Xeon Processor E5-2620 (Six Core HT, 2.0 GHz Turbo, 15 MB)	
<b>Softwares</b>	<b>Version</b>	<b>License / comments</b>
<b>Volocity</b>	6.3	Server (This computer needs to be <u>on</u> for other computer to use through server, #41710346)
<b>Metamorph</b>	7.5	Dongle #10368
<b>FIJI</b>	ImageJ 1.49	Free / updates frequently. Always use after update
<b>Zeiss Zen 2011 lite</b>	Black	Free / Viewer for .lsm, .czi from confocal
<b>Zeiss Zen 2011 lite</b>	Blue	Free / Viewer for .lsm, .czi from general Zeiss microscopes
<b>Leica LAS AF Lite</b>	3.1.0 build 8587	Free / Viewer for .lej, .lif
<b>Softworx Explorer</b>	1.3	Free / Viewer for .dv
<b>Network</b>	\\521W-130BIO21L1\Users	